520Di and 520DiN

520 series dispensing pumps




Features and benefits

- Fully automated aseptic fill/finish in biopharmaceutical and pharmaceutical drug manufacturing
- Benchtop media or reagent preparation in laboratory suites filling or inoculating bottles, vials, Petri dishes and media bags
- Five-year warranty NEMA2 and NEMA4X (IP66 and IP31) enclosure for highpressure washdown environments
- Addition of colours, dyes, flavours, vitamins, and other ingredients in food production



520Di and 520DiN performance

505L pumphead — Flow rates (ml/min)										
Tube Material				Speed	1.6mm	3.2mm	4.8mm	6.4mm	8.0mm	9.6mm
Pumpsil®, GORE® STA-PURE® Pump Tubing-Series PCS, GORE® STA-PURE® Pump Tubing - Series PFL				0.1-350rpm	0.04-150	0.23-800	0.42-1500	0.70-2500	0.90-3200	NA
Marprene®, Bioprene®				0.1-350rpm	0.04-130	0.17-500	0.38-830	0.67-2000	1.10-3200	NA
505L pumphead recommended dose sizes										
	Tube bore	1.6mm	3.2mm	4.8m	m	6.4mm	8	3.0mm	9.6mn	n
maximum speed 350rpm	Minimum dose	0.9ml	4.5ml	8.5ml		14ml	1	8ml	25ml	
	Maximum dose	15ml	75ml	150m		250ml	3	800ml	400ml	1



Dispense volume against time for the 520Di/L pump with a selection of tubes. Dispensing accuracy is indicated by the background shade




Technical specifications

	520Di and 520DiN
Accuracy	±0.5 %
Ingress rating	IP31, IP66, NEMA 2, NEMA 4X
Operating temperature range	5 to 40 °C
Operating temperature range	40 to 104 °F
Speed control ratio	2200:1
Humidity	(Condensing) 10-100 % RH, (Non-condensing) 80 % up to 31 °C (88 °F) decreasing linearly to 50 % at 40 °C (104 °F)
Noise	<70dBb(A) at 1m
Standards	CE, cETLus
Power supply	110/230 V, 50/60 Hz, Single-phase
Weight	9.82 - 10.7 kg
Weight	21.63 - 23.56 lbs
Compatible tubing size	1.6 mm ID, 3.2 mm ID, 4.8 mm ID, 6.4 mm ID, 8.0 mm ID, 9.6 mm ID

Materials of construction

	520Di and 520DiN
Bearings	Stainless steel
Casework	Pressure die-cast aluminium LM24
Keypad/HMI	Polyester
Pumphead roller assembly	Stainless steel 316
Seals	PTFE

520Di and 520DiN dimensions



Control options

Pump drives	520Di	520DiN					
Enclosure							
Ingress protection rating	IP31	IP66					
Manual control							
Digital speed control with full calibration to display flow	✓	✓					
Choice of flowrate or speed display	✓	✓					
Numerical keypad for entry of speed	✓	✓					
Remote control							
Run/Stop, direction change, auto/manual toggle, leak detector input (via contact closure or 5V TTL or 24V industrial logic)	✓	✓					
Remote operation of Dosestart with instantaneous switch	✓	✓					
Four programmable digital status outputs through 24V, 30W relays		✓					
Digital communication							
Digital network serial control	RS232	RS485					
Security							
Security code to protect set-up, programs and dose parameters	✓	✓					

Product codes

Pump and pumphead product codes							
Description	Partcode	Description	Partcode				
505L Low pulse twin channel pumphead	053.4001.000	520Di (IP31) drive only	056.7171.000*				
505LX Low pulse twin channel extension pumphead	053.4011.000	520DiN (IP66) drive only	056.717N.000**				
505LG Low pulse twin channel pumphead for STA-PURE® tubing	053.4010.00G	505LTC Tube clamp set for 505L pumphead	059.4001.000				

^{*} Supplied as 2.4m mains lead without plug. To specify American, European, United Kingdom, Australian, Argentinian or Swiss 2.8m mains lead please specify A, E, U, K, R or C respectively. Replace the last zero with the mains lead code. For example - 050.7131.10A is a 520S/R pump with American mains lead. ** Supplied as 2.4m mains cable without a plug. To specify American mains lead please replace the last zero with A.

Tube ordering codes								
Double Y tubing e	Double Y tubing elements for 520DiL dispenser							
Bore (mm)	Bioprene TL	Marprene TL	Pumpsil	Pumpsil-D	STA-PURE® PCS	STA-PURE® PFL		
1.6	933.E016.K24	902.E016.K24	913.AE16.K24	913.DE16.K24	960.E016.K24	965.E016.K24		
3.2	933.E032.K24	902.E032.K24	913.AE32.K24	913.DE32.K24	960.E032.K24	965.E032.K24		
4.8	933.E048.K24	902.E048.K24	913.AE48.K24	913.DE48.K24	960.E048.K24	965.E048.K24		
6.4	933.E064.K24	902.E064.K24	913.AE64.K24	913.DE64.K24	960.E064.K24	965.E064.K24		
8	933.E080.K24	902.E080.K24	913.AE80.K24	913.DE80.K24	960.E080.K24	965.E080.K24		
9.6	933.E096.K24	902.E096.K24	913.AE96.K24	913.DE96.K24				

Disclaimer: The information in this document is believed to be correct, but Watson-Marlow Flexicon A/S accepts no liability for any errors it contains, and reserves the right to alter specifications without notice. WARNING: These products are not designed for use in, and should not be used for, patient-connected applications. asepticsu and Accusil are registered trademarks.

wmfts.com/global



26 July 2023